

Reverse flow cycles at higher flow rates with the GORE® Protein Capture Device with Protein A



Objective

This technical note demonstrates that the GORE® Protein Capture Device (GPCD) PROA201 is capable of reverse flow cycles up to 30srt flow rates to facilitate air expulsion without significant impact to column operating pressure, flow characteristics and binding capacity.

Background

The primary design of the PROA201 GPCD is for the system to be operated in the forward flow direction (Down-Flow). However, the PROA201 GPCD is capable of limited reverse flow operation and is routinely run at 58 mL/min (60 seconds residence time) reverse flow rate as part of initial device hook-up to a liquid chromatography (LC) system.

During the initial attachment onto a LC system such as a Cytiva AKTA Pilot 600, some air may be introduced into the inlet of the device. If this occurs, then operating the system in the reverse flow direction (Up-Flow) may be performed to expel the introduced air from the system. In this situation it is desirable to flow at 80 to 100 mL/min to force the small air bubbles through the various valves and analytic instruments situated in the downstream flow path as well as flush the integrated air trap (IAT) of any bubbles for the PROA202 and PROA203 devices.

Materials and Equipment

- Cytiva AKTA Pilot 600 Liquid Chromatography System
- Four (4) PROA201 GORE Protein Capture Devices
- Chemicals outlined in Table 1
- Human Polyclonal Antibody (pIgG) with a titer of 3.0 g/L

Table 1. Chemicals utilized

| Buffer/Solution | Composition |
|---------------------------------|--------------------------------|
| Phosphate Buffered Saline (PBS) | 50 mM Phosphate 150 mM NaCl |
| Citrate | 100 mM Citrate |
| H ₂ O | |
| Sodium hydroxide (NaOH) | 0.2 M NaOH |
| 20% Ethanol (EtOH) | 20% EtOH/80% water (v/v) |

Methods

Delta column pressure

Delta column pressure was measured at 348 mL/min (10 seconds residence time) on a control device to assess initial column pressure drop. The same measurement was taken after reverse flow testing to assess any impact to column pressure performance.

Transition Test

The transition volume from PBS to water at 348 mL/min (10 SRT) was measured on three replicate devices to assess flow front performance and device integrity associated with reverse flow testing. The transition volume testing was performed per Table 2. A conductivity threshold of 0.15 mS/cm was used to determine the final volume reported.

Table 2. Transition Analysis Protocol

| Method Step | Solution/Buffer | Column Volumes (CV's) | Seconds Residence time (srt) |
|-----------------------|---------------------------------|-----------------------|------------------------------|
| Equilibration Wash I | Phosphate Buffered Saline (PBS) | 10 | 10srt |
| H ₂ O Wash | H ₂ O | 10 | 10srt |
| Equilibration Wash II | Phosphate Buffered Saline (PBS) | 10 | 10srt |

Dynamic Binding Capacity to 10% Breakthrough (DBC_{10%})

Dynamic Binding Capacity was measured at 20 SRT on the three replicate devices to assess impact to binding characteristics associated with reverse flow testing. The DBC method is outlined in Table 3 below.

Table 3. Dynamic binding capacity method

| Method Step | Solution/Buffer | Column Volumes (CV's) | Seconds Residence time (srt) |
|-----------------------|---------------------------------|-----------------------|------------------------------|
| Equilibration Wash I | Phosphate Buffered Saline (PBS) | 4 | 10srt |
| Sample Application | pIgG 2.76 g/l | As Needed | 20srt |
| Equilibration Wash II | Phosphate Buffered Saline (PBS) | 4 | 10srt |
| Elution | 0.11 mM Citrate pH-3.0 | 5 | 10srt |
| CIP Strip | 0.2M NaOH | 1 | 10srt |
| Re-Equilibration | Phosphate Buffered Saline (PBS) | ~ 8 | 10srt |

Reverse Flow Cycling

The GPCD PROA201 reverse flow capability was assessed by performing Reverse Flow, Transition test and DBC_{10%} evaluation testing sequence as described in Table 4. The test sequence was repeated for three (3) PROA201 GPCD's. Since the 58 mL/min reverse flow rate is used routinely as part of device preparation, the transition volume and DBC data represent a starting baseline. The 116 mL/min reverse flow rate is the test condition representative of the 80-100 mL/min flow rates intended for use. Delta column pressures were also measured on the devices and compared against the control.

Table 4. Reverse flow cycling progression

| Testing Sequence | Test Type |
|----------------------------|---------------------------|
| Test 1 58 mL/min (60 srt) | Reverse Flow |
| Post-Test 1 | 10 srt Transition |
| Post-Test 1 | 20 srt DBC _{10%} |
| Test 2 116 mL/min (30 srt) | Reverse Flow |
| Post-Test 2 | 10 srt Transition |
| Post-Test 2 | 20 srt DBC _{10%} |
| Post-Test 2 | Delta Column Pressure |

The GPCD must maintain ability to operate with delta column pressure (dP) ≤ 0.4 MPa. In addition, the transition volume should be ≤ 5 CV and DBC_{10%} within 10% pre and post testing. In the reverse flow, the device should not leak or burst.

Results

Delta Column Pressure Reverse Flow Results

The reverse flow delta column pressures (dP) are provided in Table 5 for reference. All are well below 0.4MPa.

Table 5. Results of reverse flow dP measurements across devices 60 srt & 30 srt

| Device | 60 srt RF dP (MPa) | 30 srt RF dP (MPa) |
|------------|--------------------|--------------------|
| RD62219-02 | 0.02 | 0.04 |
| RD62219-03 | 0.02 | 0.04 |
| RD62219-04 | 0.02 | 0.04 |

Transition Results


Transition volume and final column dP data for the three test devices are shown in Table 6 per the Table 2 test progression. The data suggest no significant impact of the 116 mL/min reverse flow testing to either transition volume or column dP as referenced by control data. All transitions are < 5 CV's.

Table 6. Results of reverse flow dP measurements across devices 60 srt & 30 srt

| Device | Post 60 srt RF TV (CV) | Post 30 srt RF TV (CV) | Equilibration |
|------------|---------------------------|---------------------------|-----------------------------------|
| | | | 10SRT Forward Flow dP (MPa) |
| RD62219-02 | 3.95 | 3.99 | 0.14 |
| RD62219-03 | 3.86 | 3.95 | 0.15 |
| RD62219-04 | 3.96 | 4.00 | 0.15 |
| RD62219-05 | Control | Control | 0.17 |

Figure 1 shows the 10 SRT transitions curves for the control device, with no reverse flow, and the three devices after both the reverse flow of 60 srt and 30 srt with no deterioration of performance noted by the presence of shoulders or shifts in the curves.

Figure 1. Chromatogram overlay of all 6 transitions curves



DBC Results

The DBC_{10%} at 20srt is within 10% when comparing the initial DBC_{10%} post 60srt reverse flow versus the DBC_{10%} post 30srt reverse flow across three devices with individual values shown in Table 7 and results relative to acceptance criteria are shown in Table 8.

Table 7. Results of DBC_{10%} Testing at 20srt across devices post reverse flow at 60 srt & 30 srt

| Device | Post 60 srt RF DBC _{10%} (g/L) | Post 30 srt RF DBC _{10%} (g/L) | DBC _{10%} Absolute Percent Difference (%) |
|------------|---|---|---|
| RD62219-02 | 37.2 | 38.0 | 2% |
| RD62219-03 | 38.7 | 38.2 | 1% |
| RD62219-04 | 37.7 | 37.4 | 1% |

Table 8. Acceptance Criteria and Summary of Results

| Acceptance Criteria | Observations | Result |
|--|---|--------|
| Transition volume ≤ 5 CV | All transition volumes ≤ 4.0 CV | Pass |
| No leaks or bursting of device | No leaks or failures | Pass |
| DBC _{10%} within 10% pre and post testing | All DBC _{10%} were within ± 2% | Pass |
| Delta column pressure (dP) ≤ 0.4 MPa | All devices < 0.4 MPa at 10srt | Pass |

Conclusions

All devices operated at delta column pressure (dP) < 0.4 MPa at 10srt. No leaks or failures were observed in any devices. All transitions were ≤ 4.0 CV and the DBC10% were within $\pm 2\%$ from initial test. The data indicate that the GPCD PROA201 is capable of operating in reverse flow up to 30 srt without loss of performance with respect to column pressure, transition volume or dynamic binding capacity. When discussing removal of air with customers, one can recommend using reverse flow at rates up to 30srt.

Gore PharmBIO Products

Our technologies, capabilities, and competencies in fluoropolymer science are focused on satisfying the evolving product, regulatory, and quality needs of pharmaceutical and bioprocessing customers, and medical device manufacturers. GORE® Protein Capture Devices with Protein A, like all products in the Gore PharmBIO Products portfolio, are tested and manufactured under stringent quality systems. These high-performance products provide creative solutions to our customers' design, manufacturing, and performance-in-use needs

NOT INTENDED FOR USE in medical device or food contact applications or with radiation sterilization. GORE Protein Capture Devices are intended for research use only and should not be used for clinical or diagnostic procedures.

All technical information and advice given here is based on our previous experiences and/or test results. We give this information to the best of our knowledge, but assume no legal responsibility. Customers are asked to check the suitability and usability of our products in the specific applications, since the performance of the product can only be judged when all necessary operating data is available. Gore's terms and conditions of sales apply to the purchase and sale of the product.

GORE, *Together, improving life*, and designs are trademarks of W. L. Gore & Associates. © 2022 W. L. Gore & Associates, Inc.

Americas | W. L. Gore & Associates, Inc.
402 Vieve's Way • Elkton, MD 21921 • USA
Phone: +1 410 506 1715 • Toll-free (US): 1 800 294 4673
Email: pharmbio@wlgore.com

Europe | W. L. Gore & Associates, GmbH
Wernher-von-Braun-Strasse 18 • 85640 Putzbrunn, Germany
Phone: +49 89 4612 3456 • Toll free: 0 800 4612 3456
Email: pharmbio_eu@wlgore.com



gore.com/pharmbio